


A male golfer is captured in the middle of a golf swing. He is wearing a black polo shirt, grey trousers, a black belt, and a black baseball cap. He is wearing a white golf glove on his left hand. A large, tan, adhesive bandage is applied to his right forearm. The background is a bright blue wall with two large, textured rectangular panels, one pink and one blue. The overall scene is brightly lit.

*Professional Healthcare
Products*



Creating Innovative Therapies



Content

- Original K-Tape* **p4**
- K-Tape My Skin* **p10**
- K-Tape Team* **p14**
- K-Tape Pure* **p18**
- K-Tape Cut* **p22**
- Crosstape* **p26**
- Accessories* **p30**
- Books* **p32**
- Biviax Worldwide* **p34**



The Original K-Tape created by Kumbrink is recognized worldwide for its optimum quality and performance.

K-Tape®... Always a Step Ahead

K-Tape is an innovative product, recognized worldwide for its optimum quality and performance. The Kumbrink developed original K-Tape combines the best natural, raw materials with the most advanced manufacturing process available.

4

- Elastic cotton tape for professional use
- Water resistant
- Sweat-permeable & breathable
- Medication free
- No active ingredients
- Physiobond acrylic glue
- 96% cotton, 4% spandex

Healthcare Professional Discount

*As professionals you are entitled to special pricing.
Visit www.k-tape.com and register.*



High-End Quality

The Cotton

Only high-quality cotton from one supplier is used for K-Tape. This quality guarantees a fabric structure with homogeneous thread density and virtually knot-free manufacture.

The Elastic Thread

If a thread loses its elasticity, the therapeutic effect is lost. This is why we use the best quality elastane thread for K-Tape. This ensures consistent elasticity levels throughout the entire period of use.

The Dyes

You need to be certain that dyes and adhesive do not harbor risks to the skin or health when applying a tape adhering to the skin for up to 7 days. We use high-quality dyes for K-Tape, supplied by Dystar. The reason: Dystar products possess the *OEKO-TEX Standard 100* label and are safety approved for use in baby clothing. K-Tape is suitable for a variety of applications and a wide range of users.

5



CERTIFIED MEDICAL PRODUCTS


MDR 2017/745

The Physiobond[®] Adhesive

Only *K-Tape* is equipped with professional grade Physiobond adhesive.

Physiobond is a registered trademark in Europe and the USA and is used exclusively for *K-Tape*. We recognize the quality of an adhesive depends not only on the choice of raw materials, but also in the manufacturing and finishing processes.

The objective during production is to create long polymer chains enabling the adhesive to anchor effectively onto the skin. This is followed by a complex temperature-time procedure to remove skin-irritating residual monomers. It is only then, that you can create a first-class tape that establishes a new standard.



The Backing Paper

Typically, backing paper is only there to protect the adhesive coating. However, with *K-Tape*, the backing paper has additional practical functions as well.

Paper Fibers

The correct orientation of paper fibers makes the backing paper more stable and *K-Tape* pre-cuts deform significantly less than those with conventional backing paper.

Prestretch

It is necessary to prestretch the cotton fabric during production to ensure crease-free manufacturing. It's an important production procedure we use for quality assurance.

K-Tape is applied to the backing paper with 10% prestretch, automatically giving it the ideal amount of prestretch for muscle technique taping. By using the correct technique to apply the tape, it is possible to transfer the required amount of prestretch to the skin. This enables *K-Tape* to always be pulled off the backing paper with the same amount of tension. Its handling always feels the same... every time you use it.

The Cut

Quite simply, it is the cut that ultimately demonstrates the outstanding quality of *K-Tape*.

A clean, straight cut not only looks professional, it also guarantees better adhesion and duration of wear. If the threads are torn by an unclear cut or if a skewed position causes the outer threads to be frequently cut, the tape unravels at a faster rate and loses its adhesiveness. We therefore take time in production, following a meticulous approach and using high-quality tools, to produce the best cut *K-Tape* possible.



Benefits for Therapist and Physician

Patients receive a high-quality product, recognized worldwide for its performance quality. Skin irritation is nearly non-existent, which means no patient complaints. Its natural materials contain no harmful substances whatsoever.

K-Tape may potentially adhere to the skin for one week. Many inexpensive tape substitutes adhere for just 2 to 3 days, meaning that they have to be changed frequently. Since patients usually have to pay for all treatments themselves, frequent renewal means higher costs for both patient and therapist.



The Original K-Tape created by Kumbrink is recognized worldwide for its optimum quality and performance. Often copied but never matched.

Single Roll (5m)

Article # 100111	1	Red
Article # 100112	2	Blue
Article # 100114	3	Black

Box of 4 Rolls (20m)

Article # 100141	4	Red
Article # 100142	5	Blue
Article # 100144	6	Black
Article # 100145	7	Mixed Colors

XXL Single Roll (22m)

Article # 100161	8	Red
Article # 100162	9	Blue
Article # 100164	10	Black





10

Invisible Natural Colors

Kitape[®]

My Skin

11

KTape[®] My Skin

K-Tape My Skin is a revolutionary skin color-matching system.

Box of 5 Rolls (25m)

Article # 100151 1 Mixed Colors

Single Roll (5m)

Article # 100115 2 Beige
 Article # 100116 3 Light Brown
 Article # 100117 4 Medium Brown
 Article # 100118 5 Dark Brown
 Article # 100119 6 Very Dark Brown

Box of 4 Rolls (20m)

Article # 100146 7 Beige
 Article # 100147 8 Light Brown
 Article # 100148 9 Medium Brown
 Article # 100149 10 Dark Brown
 Article # 100150 11 Very Dark Brown

XXL Single Roll (22m)

Article # 100176 12 Beige
 Article # 100177 13 Light Brown
 Article # 100178 14 Medium Brown
 Article # 100179 15 Dark Brown
 Article # 100180 16 Very Dark Brown



2



3



4



5



6



7



8



9



10



11



12



13



14



15



16



KITape[®] Team

Show Your Colors



Classic Sports Team Colors

Single Colors - Box of 2 Rolls (10m)

- Article # 100830 1 Sport Blue Rolls, Sport Blue Box
- Article # 100825 2 Sport Red Rolls, Sport Red Box
- Article # 100826 3 Orange Rolls, Orange Box
- Article # 100827 4 Yellow Rolls, Yellow Box
- Article # 100831 5 Green Rolls, Green Box
- Article # 100832 6 Blue Rolls, Blue Box
- Article # 100833 7 Red Rolls, Red Box
- Article # 100828 8 White Rolls, Grey Box
- Article # 100829 9 Black Rolls, Black Box

Team Color Combinations - Box of 2 Rolls (10m)

- Article # 100801 10 Sport Red & Sport Blue Rolls, Sport Blue Box
- Article # 100806 11 Black & Sport Red Rolls, Sport Red Box
- Article # 100807 12 Sport Blue & Orange Rolls, Orange Box
- Article # 100808 13 White & Orange Rolls, Orange Box
- Article # 100810 14 Yellow & Green Rolls, Green Box
- Article # 100805 15 White & Sport Blue Rolls, Sport Blue Box
- Article # 100809 16 Black & Yellow Rolls, Yellow Box
- Article # 100804 17 White & Sport Red Rolls, Sport Red Box
- Article # 100812 18 Yellow & Sport Blue Rolls, Sport Blue Box
- Article # 100802 19 White & Black Rolls, Black Box
- Article # 100803 20 White & Green Rolls, Green Box
- Article # 100811 21 Sport Blue & Blue Rolls, Blue Box
- Article # 100819 22 Sport Red & Black Rolls, Yellow Box
- Article # 100816 23 Sport Blue & Sport Red Rolls, Sport Red Box
- Article # 100821 24 Green & Sport Blue Rolls, Sport Blue Box
- Article # 100817 25 Sport Red & Yellow Rolls, Black Box
- Article # 100814 26 Yellow & Black Rolls, Black Box
- Article # 100815 27 Orange & Black Rolls, Black Box
- Article # 100813 28 Yellow & Black Rolls, Black & Grey Box
- Article # 100818 29 Sport Blue & Yellow Rolls, Yellow Box
- Article # 100820 30 Orange & Blue Rolls, Blue Box
- Article # 100822 31 Red & Blue Rolls, Grey Box
- Article # 100823 32 White & Blue Rolls, Blue Box
- Article # 100824 33 White & Red Rolls, Red Box

Customize Your K-Tape Team Box

If you do not find your team color combination in our preselected color combinations, you can customize your own K-Tape Team box by selecting 1 box color and 2 roll colors. This offer is available exclusively from biviix.



*Designed for people who
demand the purest, color-free,
100% natural cotton.*

KTape[®] Pure

Natural Color-free Cotton



K-Tape® Pure

Since its introduction, K-Tape Pure has gained considerable attention and success, as many of our customers tried the innovative new oxygen-bleached tape right away. Their reaction and approval were immediate and we received extraordinary positive feedback. They discovered its three new distinct benefits, which are a major first in the industry.

Single Roll (5m)

Article # 100120 1 Natural Cotton

Box of 4 Rolls (20m)

Article # 100152 2 Natural Cotton

XXL Single Roll (22m)

Article # 100181 3 Natural Cotton



Softer Touch

K-Tape Pure feels softer than our well-known original K-Tape. The reason? Its woven pure cotton fabric is not color-dyed. Made from tiny cotton fibers spun into a thread, the resulting tiny spaces its weave creates gives it natural elasticity and a noticeable softer feel. If the thread were dyed, those spaces would fill up with color and the fabric would lose some of its natural flexibility and feel firmer. By being 100% color-free, K-Tape Pure retains its natural flexibility and is softer to the touch. It's a significant difference that no other tape in the industry can deliver.

Stronger Adhesion

K-Tape Pure adheres even better, yet still feels light on the skin. It's made with our proprietary Physiobond adhesive, like all our K-Tape products. However, its difference is in its undyed cotton fabric. This results in an improved adhesive quality and establishes an optimal connection between the adhesive and cotton. The open-pore structure of the undyed cotton in K-Tape Pure, provides a key advantage with its unique innovative weave. The adhesive is able to penetrate much deeper into the open spaces and anchor there much better, compared with color-dyed tapes.

Longer Connection

Patients report high body-tape feedback that describes the lasting effectiveness of K-Tape Pure, even after several days. Typically, the mechanical influences from exercise, normal physical activity and showering impact tape adhesion performance. This can cause the adhesive to partially detach from the dyed-smooth thread cotton surface. However, since adhesive can better penetrate K-Tape Pure's open-space, pure cotton weave, it interlocks more firmly. The better the adhesive anchors in the cotton, the longer the skin connection lasts. This results in a longer-lasting performance with K-Tape Pure.

Ktape[®] Cut

Precision Cut. Easy to Apply.

KTape[®] Cut

K-Tape Cut was created for professional healthcare workers and their patients for easy applications.

K-Tape Cut15 (15 x 5cm) 30 Pre-Cuts

Article # 100611	1	Red
Article # 100612	2	Blue
Article # 100613	3	Beige
Article # 100614	4	Black

K-Tape Cut20 (20 x 5cm) 22 Pre-Cuts


Article # 100621	5	Red
Article # 100622	6	Blue
Article # 100623	7	Beige
Article # 100624	8	Black

K-Tape Cut30 (30 x 5cm) 15 Pre-Cuts


Article # 100631	9	Red
Article # 100632	10	Blue
Article # 100633	11	Beige
Article # 100634	12	Black



Common Applications with K-Tape Cut15




K-Tape Cut20



Shoulder & Neck

K-Tape Cut30 & Cut15




Ankle

CROSSTAPE[®]


Patches

CROSSTAPE®

The Crosstape patch is used to stimulate selected acupuncture points, enabling an acupuncture-like treatment without the use of needles. It can be used on trigger points, acupuncture points and other areas.



The 4 Crosstape patches illustrated above are shown actual size.



Small			
Article # 200100	1	400 Patches on 20 sheets	
Medium			
Article # 200101	2	180 Patches on 20 sheets	
Large			
Article # 200102	3	120 Patches on 20 sheets	
XLarge			
Article # 200103	4	40 Patches on 20 sheets	
Mix			
Article # 200224	5	Box of 6 Sheets (55 Patches)	



Common Crosstape Applications

Size M is used most of the time. Size S is ideal for facial applications, acupuncture points on fingertips, tips of toes and ears. Sizes L and XL can be used on larger areas.



For the temples, use Size M or S. For the hand, use Size M.

For the shoulder, use Size M, L or XL. For the hand use Size M.



For the elbow, use Size M, L. For the hand use Size M.



For the trigger points, use Size M or L. For the hand, use Size M.

For Stomach-Application acupuncture points, use Size M.

Professional Accessories. Professional Results.



Using professional tools, yields optimum results. Our exclusive K-Scissors have Teflon-coated blades to prevent the penetration of the acrylic adhesive into the cutting edges and to retain long-lasting sharpness. Our PreKGel® delivers superior skin preparation when applying K-Tape®.



K-Scissors K210

Article # 300101 1 Right-handers

K-Scissors K160n

Article # 300102 2 Bandage scissors

K-Scissors K130

Article # 300104 3 Mini bandage scissors with retractable clip

K-Scissors K210L

Article # 300103 4 Left-handers



PreKGel

Superior skin preparation effectively cleanses, sterilizes and stops sweat to improve adhesion, making K-Tape stick better and last longer during sports activities. Gentle, soothing and formulated for sensitive skin.

Article # 600300 5 Gel


Best-Selling Guide



Birgit Kumbrink
Author, physiotherapist,
founder and medical/
scientific director of the
international K-Taping
Academy

Birgit Kumbrink's best-selling guide to K-Taping® has become an industry standard. These richly-illustrated guide features highly successful therapeutic approaches to treating orthopedic, traumatological and many other conditions. It serves as a perfect companion to training courses and

as an important reference source and refresher for healthcare professionals, who have completed their training.



Your Go-To Reference Sources



K-Taping Illustrated Guide

Article # 603103 1 By Birgit Kumbrink (3rd edition)

K-Taping Guide illustré

Article # 604103 2 De Birgit Kumbrink (3^e édition)

K-Taping Ein Praxishandbuch

Article # 600103 3 Von Birgit Kumbrink (3. Auflage)

K-Taping in Pediatrics

Article # 603104 4 By Birgit Kumbrink (2nd edition)

K-Taping en pédiatrie

Article # 604104 5 De Birgit Kumbrink (1^{re} édition)

K-Taping bei Kindern

Article # 600104 6 Von Birgit Kumbrink (1. Auflage)

K-Taping in der Lymphologie

Article # 600105 7 Von Birgit Kumbrink (1. Auflage)

Biviax is recognized worldwide for developing its Original Kumbrink Tape Technology and creating innovative, therapeutic products of

optimum quality and performance. Biviax is active around the globe, offering superior customer service on all 5 continents.



Worldwide

Environment and Sustainability

Our environmental concerns for waste and conservation have been carefully managed through the continuing efforts of biviax, which has received the environmental certificate from Germany's Green Dot Dual System (Duales System Grüner Punkt) for *K-Tape* and *Crosstape*. The prevention of production and packaging waste, and thus the conservation of resources and reduction of greenhouse gas emissions, has continued to be a strategic objective at biviax.

Quality Control

Biviax is ISO 9001 certified for quality management. This ensures customers receive consistent, good quality products and services. Our manufacturing process meets standards set for quality and performance. This also applies to raw materials used, as the SGS test institute stringently tests our tape. Test results for *K-Tape* indicated there are no additional, unnecessary component in its fabrication.

K-Taping® International Academy


Professional training is essential in order to effectively use *K-Tape*. Our certified training program courses, created by Academy founder and medical director Birgit Kumbrink, are offered worldwide. View *K-Taping* courses available or register for a course in your region at www-k-taping.com

Want to Become a Distributor?

There are many benefits to becoming an authorized distributor of *K-Tape* products. Not only is there exclusivity, but also the opportunity to grow your business and image through association with the biviax global brand. For detailed information, please contact us directly at info@biviax.com

New

Scissors® K130



Mini Bandage Scissors

Retractable Clip



Creating Innovative Therapies

Phone +49 (231) 9767-301 • Fax +49 (231) 9767-302 • info@biviax.com • www.biviax.com • www.k-tape.com
biviax GmbH & Co. KG, Hagener Strasse 428, 44229 Dortmund, Germany

© K-Tape, K-Tape My Skin, K-Taping, Physiobond, Crosstape, Crosspuncture and PreKGel are registered trademarks in Germany, Europe, USA, Canada and many more countries. All rights reserved.